# 383 nm Laser system (T-MOT)

The 383 nm laser system consists of 767 nm lasers which are then frequency doubled to achieve 383 nm light.

### **External Cavity Diode Laser (ECDL)**

For 767 nm lasers, ECDL in Littrow configuration is used. Typically we used the laser diodes from Eagleyard Photonics: EYP-RWE-0790-02000-1500-SOT02-0000

Recently, Eagleyard has replaced these with new laser diodes: EYP-RWE-0760-02010-1500-SOT12-0000

#### TA

Output Power: 1.5 WInput Current: 2 A

• Injection Power: 32 mW

Power behind 30dB Isolator: 1.05 W

• Originally this TA was used: EYP-TPA-0765-01500-3006-CMT03-0000. Is this still true?

#### **Fiber**

• PMC-780-5,0-NA012-3-APC-200-P

• Incoupling: 67%

• Power behind fiber: 700 mW

## Frequenzy doubling

### **LBO-Crystal**

• Lenght: 15 mm

AR coating

#### Resonator

• Ring resonator (double Z configuration)

• Length: 280mm

• Curvature of mirrors: 50 mm (S3 and S4)

· Distance of mirrors: 64 mm


Waist: 30μm (crystal), 130μm (long arm)

• Transmission: TS3 = 0.049%, T1 = 1.2 %

• Conversion efficiency ENL = 6.1\*10^(-5)/W


 $\$E \{NL\} = \kappa L c k 1 h m (B,\kappa) $$$ 

- Linear losses: eL = 0.85(0.15) %
- Finesse: F = 270




# **Stabilisation**

- PDH-Method
- Error signal at about 20 MHz




# Stabilisation: Laser 4 - Potassium

Dopplerfree Saturation Spectroscopy on D1 line of Potassium:



- 1st the two beams with similar intensity are generated by the beam sampler and are send trough the glas cell
- one part of the initial beam is going trough the sampler to a double pass aom and then in the other direction trough the glas cell, crossing only one beam. Important: the beam coming from the aom must have a much higher intensity!
- the first two beams, coming from the sampler trough the cell are than monitored by the PD. The signal is substracted
- in or case e use the aom to create the peak on one side of the mainpeak. This is important for the lock-in scheme



From:

https://iqwiki.iqo.uni-hannover.de/ - IQwiki

Permanent link:

https://iqwiki.iqo.uni-hannover.de/doku.php?id=groups:mg:experiment:laser:383nmlaser

Last update: 2022/04/01 09:07

