


Notes/Informations

Frequency-Comb Hannover

Channel	Measured	Frequency	Sign	Interval	Max outliers (Alarm)	F(x)
CH1	IQ-Resonator (R1) vs. Comb	29.347 MHz	-	± 1 kHz		-fx
CH2	CEO	20 MHz	-	± 2 Hz		-fceo
CH3	Reprate	6.250001 MHz		± 0.002 Hz		frep
CH4	Link vs. Comb	22.432 MHz	-	± 300 Hz		-fy
CH5	Virtualer Beat (Link vs. IQ-Laser)	21,27143 MHz		± 20 Hz		fv
CH6	Manypoint-AOM	36 MHz		± 5 Hz		
CH7	Stell-AOM vom IQ-Laser	35.02865 MHz		± 100 kHz		
CH8	Rubidiumuhr referenziert auf GPS	10 MHz		± 0.002 Hz		

- Modenverhältnis: 1944000/3275292
- Repetitionsrate: 100,0000226 MHz



- AOM (+77.111 MHz) used for atom-lock + fiber stabilization(5m)
- AOM (+80.111 MHz) used for fiber stabilization(30m)
- Double-pass AOM (-2*102,XX MHZ) used for spectroscopy (frequency shift), feedback-loop!!!

- Channel 7 = 138 MHz (Synthesizer, referenced by 10 MHz GPS) - Double-pass AOM

Virtuel-Beat

- $f_v = (f_y + f_{ceo}) - (f_x + 2*f_{ceo}) * (1944000/3275292) + (2*10 \text{ MHz})$

Evaluation (PTB)

- $f_x = \text{col}(\text{C}) \#1$, negatives f_x
- $f_{ceo} = \text{col}(d) \#2$ (negatives f_{ceo})
- $f_{rep} = \text{col}(e) \#3$, $f_{rep}/16$

Test analysis

- $\nu_{opt}(\text{Mg}) - \nu_{\text{(Mg2015)}}[\text{Hz}] = \nu_{opt}(\text{Mg}) - 655058646691000$
- Atom frequency: $\nu_{opt}(\text{Mg}) = (3275292*\text{col}(e)*16 - 2*\text{col}(d) - \text{col}(\text{C}) + 80,111\text{e}6)*2 + 77,111\text{e}6$
- Laserfrequency: $f_{laser} = 3275292*f_{rep}*16 - 2*f_{ceo} - f_x$

→ $\nu_{opt}(\text{Mg}) - \nu_{\text{(Mg2015)}}[\text{Hz}] = (f_{laser} + 80,111\text{e}6)*2 + 77,111\text{e}6 - 655058646691000$

Frequency-Comb Braunschweig

From:

<https://iqwiki.iqo.uni-hannover.de/> - IQwiki

Permanent link:

<https://iqwiki.iqo.uni-hannover.de/doku.php?id=groups:mg:notes&rev=1490707841>

Last update: **2017/03/28 13:30**

