


# Ring-down Box

## Application of the ring-down box

You can measure the finesse of a resonator with the ring-down box **without** locking the resonator!

### Advantages

1. super fast (see characterization)
2. without locking the resonator



### Measurement Method

### Evaluation

Use this formula to analyse the datapoints:  
 $y(x) = A \cdot \exp(-x/t) + y_0$  with  $x := \text{time}$ ,  $t := \text{delay time}$

Afterwards calculate the finesse with  $F = 2\pi \cdot t \cdot \text{FSR}$  with  $\text{FSR} := \text{free spectral range}$ .


### Electronic setup


- max. Input PD-voltage: +5V or -5V
- Output TTL (Name: inv. Q): 5V
- Power voltage: +/-15 or +/-18V (max.)
- Switch A: unten = Eingang invertierend, oben = deaktiviert

- Switch B: *unten = Ringdown Messung, oben = dauer-reset*

## **Explanation of the operating principle**

- Electronic drawing (PDF-file)





## Characterization

### List of components

- Inverter: LT1363
  - Invert negativ input voltage to positiv
    - Datasheet\_LT1363
  - Power supply: +/-5V
- Comparator: AD8561
  - Compare input voltage and threshold value
    - Datasheet\_AD8561
  - Power supply: +/-5V
  - High speed timing: 7ns @ 5V
- Flip-Flop: SN74HC74
  - Output of trigger
    - Datasheet\_SN74HC74
  - Power supply: +5V



- High speed timing: 15ns
- 1x XLR connector
- 1x Power supply +/-5V
- 1x single switch
- 1x double switch
- 3x BNC connectors
- 1x Potentiometer
- 1x big red button
- 2x 8-pin holder
- 1x 14-pin holder
- 3x PSK 3-pin cables
- 5x 100n capacitors
- some resistors (3x 100k, 1x 4.7k, ...)
- 1x Al-box
- 2x PCB's

## Microwave switch

- Control of RF signal to the AOM
- Mixer typ: SRA-1+ ( [Datasheet\\_SRA-1+](#) )
- 1x Resistor (100 Ohm)
- 2x diodes (typ: 1N4748)
- 1x diode (typ: Z3,9V)
- 3x BNC connectors
- 1x Al-Box



But all ground pins on ground, also the housing



From: <https://iqwiki.iqo.uni-hannover.de/> - IQwiki

Permanent link: [https://iqwiki.iqo.uni-hannover.de/doku.php?id=groups:mg:ring-down\\_box&rev=1516882927](https://iqwiki.iqo.uni-hannover.de/doku.php?id=groups:mg:ring-down_box&rev=1516882927)

Last update: 2018/01/25 12:22

