



Fiberlaser


- Push Yes



- Push OK



- If you **turn off and on** the frequency comb, then push in the order:
 1. *Reset*
 2. *Auto on*
 3. *Enable on*
- For **closing** the software:
 1. *Halt*
 2. *Freeze*
 3. Close the window



From:
<https://iqwiki.iqo.uni-hannover.de/> - IQwiki

Permanent link:
<https://iqwiki.iqo.uni-hannover.de/doku.php?id=groups:mg:fiberlaser&rev=1481624543>

Last update: **2016/12/13 10:22**

