



# List of Computer-Programs



Save a backup file from each program in AFS!

| Name                                                                                                                               | Description             | Computer | Information                                                                                |
|------------------------------------------------------------------------------------------------------------------------------------|-------------------------|----------|--------------------------------------------------------------------------------------------|
|  <a href="#">TEC Service XY</a>                   | TEC Control Software    | "        | Meerstetter Control Software for ECDL baseplate, active shieldings of resonators R1 and R2 |
|  <a href="#">Allan-Live with time graph</a>       | -                       | "        | -                                                                                          |
|  <a href="#">TracerDAQ</a>                        | RedLab Card Counter     | ResLab   | Can count Voltages; Limited Capacities / Time;                                             |
| <a href="#">RbMon</a>                                                                                                              | Rubidium clock software | "        | -                                                                                          |
|  <a href="#">Net View</a>                        | GPS receiver            | "        | Control software for GPS signal.                                                           |
|  <a href="#">Fiber Comb Control</a>             | ...                     | "        | This Software can control Repetition Rate Stepper Motor and CEO Tuner                      |
|  <a href="#">FiberLaser</a>                     | First start this        | "        | First start this, better Autotuning                                                        |
|  <a href="#">Kamm-Frequenzneu kombiniert.nb</a> | -                       | "        | -                                                                                          |
|  <a href="#">allen-combined</a>                 | -                       | "        | -                                                                                          |
|  <a href="#">Spyder</a>                         | Pyhton Compiler         | "        | Code for Einstein-Telescope                                                                |
| Pyhtonberechnungen von und mit Schlippi Folder                                                                                     | -                       | "        | -                                                                                          |
|  <a href="#">K&amp;K FXE / FXQE80</a>           | Counter                 | "        | Maximal 8 Com Ports                                                                        |
|  <a href="#">KK FXE Folder</a>                  | Folder of Counter       | "        | Counter Daten werden hier gespeichert                                                      |
| ...                                                                                                                                | ...                     | ...      | ...                                                                                        |
| ...                                                                                                                                | ...                     | ...      | ...                                                                                        |

| Name                                                                                                      | Description                  | Computer     | Information                                                                                                                                                                                              |
|-----------------------------------------------------------------------------------------------------------|------------------------------|--------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  <b>Monitor VCH-1006</b> | Control Software for H-Maser | ResLab Maser | one entry is always red! Andreas Bauch doesn't know the reason, but it's okay!                                                                                                                           |
|  <b>DDS USB</b>          | Programming the USB-DDS      | ResLab Maser | -                                                                                                                                                                                                        |
|  <b>USB Serial</b>       | Programming the USB-Serial   | ResLab Maser | -                                                                                                                                                                                                        |
| ...                                                                                                       | ...                          | ...          | ...                                                                                                                                                                                                      |
| <b>Highfinesse Wavemeter WS/6-200</b>                                                                     | Wavemeter                    | ResLab Mobil | Push the "Start Measurement" button                                                                                                                                                                      |
| Laser Frequency ResLab.vi                                                                                 | Wavemeter                    | "            | Shows Laser Frequency, also for Atom Lab                                                                                                                                                                 |
| <b>K&amp;K FXE / FXQE80</b> (currently not in use!)                                                       | Counter                      | "            | Default Settings: 8 Com Ports, Trimmer: Auto Trimmer, Measurement Method: Phase; Push buttons before measurement to erase RAM: Clear Buffer, Reset; Start Measurement via Häckchen setzen auf frequency; |
| <b>KK FXE Folder</b> (currently not in use!)                                                              | Folder of Counter            | "            | Counter Daten werden hier gespeichert                                                                                                                                                                    |
| ...                                                                                                       | ...                          | ...          | ...                                                                                                                                                                                                      |

From:  
<https://iqwiki.iqo.uni-hannover.de/> - IQwiki

Permanent link:  
<https://iqwiki.iqo.uni-hannover.de/doku.php?id=groups:mg:listprograms&rev=1476087939>

Last update: **2016/10/10 08:25**

