

Ultra-stable local oscillator (H-Beast)



- Labbook
- ToDo
- List of ideas/open questions
- Wishlist
- List of orders and status
- Paper
- List of collaboration partner and companys

Setup

Resonator

- Mirrors/Contamination
 - Substrates
- **Finesse measurement**
 - **Finesse measurement outside vacuum**
 - Finesse measurement without spacer
 - Finesse measurement with spacer
 - Birefringence measurement with spacer
 - **Finesse measurement inside vacuum**
- ULE-compensations-rings
- Spacer
- Optical Mask
- Optical Contacting
- **Calculations**
 - Calculation of Finesse outside vacuum
 - Calculation of reflectiondip / Impedance of the cavity
 - Calculation between cavity spectrum and frequency distance
 - Calculation of EOM frequency

- Calculation of linewidth

Mounting

- Mounting
- Cleaning of the mounting

Vacuum

- Vacuum tests
 - Vacuum tests
- Bake out
 - Bake out tips
 - Bake out List
 - Bake out protocol
- Vacuum Components
 - Vacuum chamber
 - CF-components: Pump, Sensor, Valves, etc.
 - Attaching of the vacuum feedthroughs
 - Seal material and vacuum suitable glue
 - Anpressring für die Fenster
 - Heat shield's
 - Components: Heat shield and vacuum suitable screws
 - Electrolytic polishing
 - Mass spectrometry
 - Determination of transport screws for the resonator
 - Attaching of the windows on the heat shields
 - Vacuum optics
 - Optical components inside the vacuum: Windows and isolator
 - Attaching of the windows on the vacuum chamber

Optics

- Incoupling breadboard
 - Incoupling breadboard
 - Components
 - Mode Matching
- Outcoupling breadboard
 - Outcoupling breadboard

Electronics

- **Electronics inside the vacuum**

- Components for the electronics inside the vacuum and accessories
- Instructions for attaching the electronic and labelling of the NTC's
- Characterization of the NTC's


- **Electronics outside the vacuum**

- PTB electronic assignment
- Components for the electronics outside the vacuum
- Labelling of the cables

Laser at PTB

- Company: NKT Photonics
- Fiber laser Koheras Basik Module
- Product type: K122-102-665
 - Wavelength: 1542.14 nm
 - Power setpoint: 40 mW
 - Power Range: 12-40 mW
 - Wavelength tuning: -499.9 ... +570 pm
 - Power stability 0.01 dB

Laser at IQ



- Todo List
- Wavelength: 1560.47 nm
- ITU channel number 21, 192.10 THz, 1560.61 nm
- [Laserdiode](#)
- [Laserdiode Mount](#)
- [Laser driver](#)
- [Pin assignment/Cables](#)
- [Temperature Controller](#)
- [Laser Housing](#)
- [Characterisation](#)

Fiber setup

- Fiber-components
 - Fiber-Photodiode (IQ Version)
 - EOM and AOM Box
- Electronic-Components
- Fiber-connections
- **Vibration isolation**


Performance / Measurements

- Thermal noise limit
- Pound-Drever-Hall (PDH)
- Residuelle Amplitudenmodulation (RAM)
- Vacuum pressure

Future Setup @ IQ

• Transfer of stability:

- Providing 10–16 Short-Term Stability of a 1.5- μ m Laser to Optical Clocks
, C. Hagemann et. al., IEEE Transactions on instrumentation and measurement, VOL. 62, NO. 6 (2013)



From:
<https://iqwiki.iqo.uni-hannover.de/> - IQwiki

Permanent link:
https://iqwiki.iqo.uni-hannover.de/doku.php?id=groups:mg:project_ptb-cavity&rev=1520497167

Last update: 2018/03/08 08:19

