




Checklist Resonator Lab (D215)


1. a) Check **IGP Controller** of resonator R1, Current < 2 mA
(should show „HV ON“ (green LED on))
b) Check both **USV for IGP**




2. a) Check both **USV's for Meerstetter TEC** of R1
(green LED's on)
b) Check **TEC** for ECDL (green LED on)



3. Check the H-Maser
(3 green LEDs on)



4. a) Check **USV** for GPS (green LED on)
b) Check **USV** for 10 MHz (green LED on)



In emergency, please call:

Steffen Sauer:

Klaus Zipfel:

Nandan Jha:

Waldemar Friesen-Piepenbrink:

0172 2068556

0160 1602392

0151 6357240

0176 30778529

- Temperature 21°C?
- See water?
- Smell burning/see smoke?
- Computer are running?