

How to flush the Coils

Preperation

Before you start to the procedure you should get everything together, that is needed.

You need

- Vinegar essence
- Destilled Water
- Boxes (I recommend at least two)
- Submersible pump
- Power supply
- Blind plugs 2x
- Spanner (19, 17 & 14)
- Some paper

Getting started

1. ...

From:
<https://iqwiki.iqo.uni-hannover.de/> - IQwiki

Permanent link:
<https://iqwiki.iqo.uni-hannover.de/doku.php?id=groups:mg:tips:coils&rev=1593603594>

Last update: **2020/07/01 11:39**

